

GST Pulldown – for 6 well plate

Transfect cells with transfection reagent of choice.

1. Wash cells with 500 μ l 1x PBS
2. Lyse cells by adding 350 μ l IP/Lysis Buffer (with 0.1% NP-40 and add PIC fresh)
3. Incubate on ice for 30 min., vortexing every 10 min.
4. Centrifuge 14,000 rpm, 4°C, 15 min.
5. Save 80 μ l lysate for total lysate and add to 80 μ l 2x Tris/tricine sample buffer
6. Prepare Glutathione Sepharose beads (20 μ l beads/ sample) by washing total amount of beads 2x with 500 μ l 1x PBS
 - a. Centrifuge 500 x g for 5 min. to pellet beads
 - b. Pipet off buffer, then add 1x PBS and proceed with washes
 - c. Centrifuge in between washes 500 x g for 5 min.
 - d. Pipet off supernatant
 - e. After final wash, resuspend beads in 1x PBS (20 μ l 1x PBS/ number of samples)
 - f. Aliquot 20 μ l washed beads into 1.5 ml microcentrifuge tubes
7. Add remaining lysate to 1.5 ml microcentrifuge tube containing 20 μ l Glutathione Sepharose beads
8. Rotate in cold room for 2 hours
9. Centrifuge 500 x g for 5 min., 4°C
10. Wash 4x with 250 μ l IP/Lysis buffer (with 0.1% NP-40 and PIC fresh)
 - a. Centrifuge in between washes 500 x g, 4°C for 5 min.
 - b. Pipet off supernatant
11. After final wash, pipet off most of the supernatant and add 20 μ l 2x Tris/tricine sample buffer to a total of approximately 40 μ l.

Lysis Buffer:

137 mM NaCl
10 mM Tris, pH 7.4
1% NP-40

2x Tricine Sample Buffer (for 50 ml):

5 ml 1M Tris, pH 6.8
12 ml Glycerol
4 g SDS
1.55 g DTT
10 mg Coomassie Blue R250
to 50 ml with dH₂O