

TOPO cloning

PCR cloning

For a 50 μ L reaction using Stratagene Pfu Turbo DNA Polymerase (cat# 600250, blunt ends):

1. Combine the following in a PCR tube, adding Taq last.

5 μ L	10x reaction buffer
2 μ L	dNTPs 5 mM each
1 μ L	Forward primer (200 ng/ μ L)
1 μ L	Reverse primer (200 ng/ μ L)
35.8 μ L	dH ₂ O
0.2	Pfu Turbo
5 μ L	template (10ng/ μ L)

2. Place samples in thermocycler using the following program:

95°C	2 min
95°C	30 sec
T _m *-5°C	30 sec
72°C	1 min/kb
72°C	10 min
4°C	hold

*T_m is calculated as $2x(\#A + \#T) + 4x(\#G + \#C)$

3. Run samples on an agarose gel to verify amplification. If other bands are present other than the expected it will be necessary to gel purify the correct fragment (QIAGEN cat#28704).

Adding A's

For gel extracted DNA:

1. Combine the following to the gel extracted fragment (for 50 μ L eluate).

5 μ L	10x Reaction buffer
2 μ L	dNTPs 5 mM each
1 μ L	Taq DNA Polymerase

2. Incubate at 72°C for 10 min.

For PCR product:

1. Combine the following to the remaining PCR product.

1 μ L	Taq DNA Polymerase
-----------	--------------------

2. Incubate at 72°C for 10 min.

TOPO TA cloning (Invitrogen cat#K2500-20)

1. Combine the following:

4 μ L	PCR product
1 μ L	Salt solution
0.5 μ L	TOPO vector

2. Incubate at RT for 10 min.
3. Transform 2 μ L into 1 vial TOP10 Oneshot competent cells.
4. Incubate transformation on ice for 30 min.
5. Heat shock transformation for 45 sec at 42°C.
6. Incubate on ice for 2 min.
7. Add 250 μ L SOC and incubate with shaking 225 rpm, 37°C for 1 hour.
8. Plate 200 μ L transformation on LB-Spectinomycin (100 μ g/mL Spec) plates.
9. Incubate overnight at 37°C.