

Mutagenesis Using FailSafe Buffers with Stratagene Kit

- 1) Use the FailSafe PCR 2X PreMixes which produced the correct size fragments during the Selection Procedure. Set up the following rxn, on ice:
 - X μ L plasmid template (30 ng)
 - X μ L primer 1 (125 ng)
 - X μ L primer 2 (125 ng)
 - 1 μ L Epicentre FailSafe PCR Enzyme Mix
 - X μ L to sterile water (to 25 μ L)Add:
 - 25 μ L FailSafe PCR 2X PreMix (A-L, use one that worked)
- 2) Run samples in the thermocycler according to the following parameters:
 - 95°C – 30s 1 cycle
 - 16 cycles: 95°C – 30s
 - 55°C – 1 min
 - 68°C – 1min/kb of plasmid length
 - 4°C hold
- 3) Add 1 μ L *DpnI* restriction enzyme (Stratagene Kit, -20°C) directly to each PCR rxn. (*DpnI* digests methylated DNA, leaving the newly synthesized mutated DNA untouched.) Mix gently. Incubate at 37°C for 1 hour.
- 4) Thaw XL-1 Blue Supercompetent cells (Stratagene Kit, -80°C) on ice. Aliquot 50 μ L of cells into pre-chilled PCR tubes.
- 5) Transfer 1 μ L of DpnI-treated DNA from each PCR reaction to an aliquot of cells. Stir with pipette tip.
- 6) Incubate on ice for 30 min.
- 7) Heat pulse the transformation rxns for 45s at 42°C (use the thermocycler). Incubate on ice for 2 min.
- 8) To culture tubes, add 0.5 mL of S.O.C. Medium (pre-warmed to 37°C). Pipette transformed cells into the media. Shake for 1 hour at 225 rpm.
- 9) Plate 200 μ L of each transformation reaction onto plates. Incubate O/N at 37°C.